

## Goat Anti-Rat IgM Antibody

### ORDERING INFORMATION

Catalog No.	Size
88601	1.0mg

**Format:** 1mg/ml in phosphate-buffered saline (PBS), 0.09% sodium azide.

### SPECIFICATION SUMMARY

**Antigen:** Rat IgM

**Host Species:** Goat

**Specificity:** Antibody is affinity-purified on rat IgM. This antibody may cross-react with IgM of other species.

### APPLICATIONS

*Immunoblotting:* use at a dilution of 1:1,000-1:30,000.

*Immunohistochemistry:* use at a dilution of 1:200-1:2,000.

*Immunocytochemistry:* use at a dilution of 1:200-1:2,000.

*ELISA:* for detection use at a dilution of 1:1,000-1:30,000; for coating use at a dilution of 1:100-1:500.

These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.

### DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

### STORAGE AND STABILITY

This antibody is stable for at least 1 year at 2-8°C.

*For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.*