Product Data Sheet Recombinant Antigen



Product name	HSV-1 Glycoprotein G	
Order no.	BA1051RVS-G	
Lot no.	Lot specific	
Source	E. coli derived recombinant protein	
Product description	The protein contains an internal fragment of HSV-1 glycoprotein G.	
Purification	By chromatography	
Presentation	The antigen is presented in buffer containing 50 mM Tris-HCl pH 8.0, 150 mM NaCl and 10 % glycerol.	
Protein concentration	Lot specific, specified in mg/mL	Determined by Bradford assay (BioRad) using BSA as standard.
Analysis	By SDS-PAGE and ELISA	
Storage	-65 °C or lower	Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific	
Expiry date	10 years from date of manufacture (under recommended storage conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.