

Product Data Sheet Native Antigen



Product name	Legionella pneumophila Antigen	
Order no.	BA106VS	
Lot no.	Lot specific	
Agent	Pool of Legionella pneumophila serovars 1A, 1B, 2, 3, 4, 5, 6 and 7	
Culture system	Solid medium	
Preparation	Washed and heat treated bacteria are presented as a suspension in phosphate buffered saline.	
Inactivation	Heat treatment. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.	
Preservatives	None	
Protein concentration	Lot specific, specified in mg/mL	Determined by Bradford assay (BioRad) using BSA as standard.
Analysis	By ELISA	
Storage	-65 °C or lower	Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific	
Expiry date	10 years from date of manufacture (under recommended storage conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area with soap and water. If irritation occurs, consult a physician.