## **Product Data Sheet Recombinant Antigen**



Product name	Parvovirus VP2
Order no.	BA122VSVP1E
Lot no.	Lot specific
Source	E. coli derived recombinant protein
Product description	Parvovirus B19, recombinant VP2 (full length, linear epitopes, VP1 aa 233-781). Molecular mass: 63 kDa
Purification	By chromatography
Presentation	The protein is presented in buffer containing 50 mM Tris/HCl pH 8.0, 500 mM NaCl and 8 M urea.
Protein concentration	Lot specific, Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By SDS-PAGE
Storage	-65 °C or lower Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from date of manufacture (under recommended storage conditions)
Manufacturer	Institut Virion\Serion GmbH, Germany

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.