Product Data Sheet Recombinant Antigen



Product name	Parvovirus VP2 VLP	
Order no.	BA122VSVP2	
Lot no.	Lot specific	
Source	Recombinant protein, expressed in SF9 insect cells	
Product description	Parvovirus B19, recombinant VP2 VLPs (VP2 virus like particles, conformational epitopes, VP1 aa 228-781)	
Purification	By chromatography	
Presentation	The antigen is presented in phosphate buffered saline at pH 7.4, containing 10% glycerol.	
Protein concentration	• •	nined by Bradford assay (Biorad) BSA as standard.
Analysis	By SDS-PAGE and ELISA	
Storage	-65 °C or lower Avoid r	epeated freeze-thaw cycles.
Date of manufacture	Lot specific	
Expiry date	10 years from date of manufacture (under recommended storage conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.