Product Data Sheet Recombinant Antigen



Product name	Mycoplasma pneumoniae P1	
Order no.	BA127VSP1	
Lot no.	Lot specific	
Source	E. coli derived recombinant antigen	
Product description	Mycoplasma pneumoniae, recombinant adhesin protein P1. Molecular mass: 37 kDa	
Purification	By chromatography	
Presentation	The protein is presented in buffer containing 50 mM Tris/HCl pH 8.0, 500 mM NaCl and 8 M urea.	
Protein concentration	Lot specific, specified in mg/mL	Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By SDS-PAGE	
Storage	-65 °C or lower	Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific	
Expiry date	10 years from date of manufacture (under recommended storage conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.