

Product Data Sheet Native Antigen



Product name	Aspergillus fumigatus Antigen	
Order no.	BA132VS	
Lot no.	Lot specific	
Agent	Aspergillus fumigatus, strain Ag-507	
Culture system	Broth culture	
Preparation	Disrupted somatic material clarified by centrifugation in phosphate buffered saline.	
Inactivation	Viable organisms are removed by filtration. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.	
Preservatives	None	
Protein concentration	Lot specific, specified in mg/mL	Determined by Bradford assay (BioRad) using BSA as standard.
Analysis	By ELISA	
Storage	-65 °C or lower	Avoid repeated freeze-thaw cycles. Sonicate before use.
Date of manufacture	Lot specific	
Expiry date	10 years from date of manufacture (under recommended storage conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.