## **Product Data Sheet Recombinant Antigen**



Product name	EBV Capsid Antigen P23
Order no.	BA1361VSR21
Lot no.	Lot specific
Source	E. coli derived recombinant protein
Product description	Epstein-Barr virus (EBV), recombinant capsid antigen P23. Molecular mass: 21 kDa
Purification	By chromatography
Presentation	The protein is presented in buffer containing 50 mM Tris/HCl pH 6.0, 1 M NaCl and 4 M urea.
Protein concentration	Lot specific, Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By SDS-PAGE and ELISA
Storage	-65 °C or lower Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from date of manufacture (under recommended storage conditions)
Manufacturer	Institut Virion\Serion GmbH, Germany

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.