Product Data Sheet Recombinant Antigen



Product name	EBV Early Antigen P54	
Order no.	BA1363P54VS	
Lot no.	Lot specific	
Source	E. coli derived recombinant protein	
Product description	The protein contains a fragment of EBV (Epstein–Barr virus) early antigen p54. Molecular mass: 45 kDa	
Purification	By chromatography	
Presentation	The antigen is presented in phosphate buffered saline at pH 7.4, containing 4 M urea.	
Protein concentration	Lot specific, specified in mg/mL	Determined by Bradford assay (Biorad) using BSA as standard.
Analysis	By SDS-PAGE and ELISA	
Storage	-65 °C or below	Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific	
Expiry date	10 years from date of manufacture (under recommended storage conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.