

Product Data Sheet Recombinant Antigen

Product name	Chlamydia trachomatis MOMP
Order no.	BA1372VSR2
Lot no.	Lot specific
Source	E. coli derived recombinant protein
Product description	Chlamydia trachomatis, fragment of MOMP (major outer membrane protein). Molecular mass: 23 kDa
Purification	By chromatography
Presentation	The protein is presented in buffer containing 50 mM Tris/HCl pH 8.0, 500 mM NaCl and 8 M urea.
Protein concentration	Lot specific, Determined by Bradford assay (Biorad) specified in mg/mL using BSA as standard.
Analysis	By SDS-PAGE and ELISA
Storage	-65 °C or lower Avoid repeated freeze-thaw cycles.
Date of manufacture	Lot specific
Expiry date	10 years from date of manufacture (under recommended storage conditions)
Manufacturer	Institut Virion\Serion GmbH, Germany

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area thoroughly with soap and water. If irritation occurs, consult a physician.