

## serion\immunologics Product Data Sheet - Native Antigen

Name of product	Parainfluenza Virus Type 3 Antigen	
Order No.	BA1263VS	
Lot No.	Lot specific	
Agent	Parainfluenza virus type 3, strain C 243	
Culture system	Vero cell culture	
Preparation	Virus is extracted and separated from host cells and cell culture supernatant, further purified by size exclusion chromatography and concentrated. The antigen is presented in Tris/NaCl buffer at pH 7.4.	
Pathogen inactivation	β-Propiolactone treatment.  No specific biological hazards are known to be associated with this product. It should be handled with appropriate safety precautions for biological materials.	
Storage	-65 °C or below. Avoid repeated freeze-thaw cycles. Sonicate before use.	
Date of manufacture	Lot specific	
Expiry	10 years from date of manufacture (when stored under recommended conditions)	
Manufacturer	Institut Virion\Serion GmbH, Germany	

For research or manufacturing use only. Not for use in diagnostic procedures.

## **Analysis:**

Test	Method	Specification
Protein concentration	Bradford assay (Bio-Rad) with BSA as reference standard	Report in mg/mL
Immunoreactivity	SERION ELISA classic Parainfluenza Virus IgG (ESR126G)	QC serum panel (multiple negative, borderline and positive sera) within defined reactivity range
Inactivation assay	Attempted pathogen recultivation under original culture conditions	No propagation observed